



PRODUCT
CATALOG

BioLine[®]

Dental Implants Series



THE ART OF SIMPLICITY

 **Notes:** _____

All rights are reserved.
no part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any form or by any means, electronic, mechanical, photocopying, recording or otherwise, without prior permission of the publisher.
warning: only a licensed dentist should use these products.

| | |
|--|---------|
| About BioLine | 4 |
| R&D departments | 5 |
| Quality system | 6 |
| Labeling & packaging | 7 |
| Bio <i>fix</i> surface | 8 |
| Case report | 9 - 10 |
| Drilling protocol | 11 |
| BioGuide | 12 |
| Bio materials | 13 |
| Internal hex implants | 14 - 27 |
| The tilt concept | 28 - 38 |
| Internal hex prosthetics | 39 - 51 |
| Conical connection dental implants | 52 - 55 |
| Conical connection prosthetic's | 56 - 59 |
| Instruments & drills | 60 - 63 |
| Surgical kits | 64 - 68 |
| Implant motor | 69 |
| Intraoral scanner | 70 |
| Digital components | 71 |

The art of simplicity, is our guiding line

At Bioline Dental Implants, we address the needs of dental implantology's, sharing our experience and technologies, utilize experts and invest in research, training and education. All this results in the art of simplicity solutions that simplify dental implantology procedures and deliver proven clinical long-term success.

We maximize our experience and innovative technologies

We Maximize our experience and technologies to ensure that our products offer the best high-end clinical solution. From the beginning of our activity in 2005, we have been focused on the development, design and manufacture of dental implants and very accurate prosthetic parts. Our state of the art manufacturing facility, which is operational 24/7, includes a dedicated team of designer's clinical consultant, cnc operators and Qa department to ensure the highest possible standards and quality of our products enable us to grant a lifetime warranty for our dental implants.

We offer high quality, simple-to-use product lines

We have been focused on state of the art simple-strong-accurate implant technology by developing implants based on two connections principals - an internal hex connection and a conical hex connection. Created on a switching -platform that fit all possible prosthetic system, they allow dentists to use a single restoration line for each implant platform. Furthermore, the development of all our products –implants, abutments and surgical tools – takes existing systems into account in order to make the procedure simple and minimize dental surgent period of learning. As a resolute of this method, our surgical tool kit is compatible with all products and includes everything from basic surgical tools to advanced instruments needed in maximum efficiency and agronomic – economic thinking. And in addition to make work processes short and simple, our products have an overall implant clinical success rate of 98.7% reported and documented according to our detailed recommended instruction for use.

We invest in research and development constantly

Our R&D teams collaborate closely with an international panel of experts who have extensive clinical and academic knowledge. We also invest in new technologies and encourage innovative thinking, clinical trials, studies. This investment of research and development ensures that our products are constantly aligned with market needs and constant improvement's, and that we are able to develop clever, easy and simple to use, best quality products.

We Invest in Training and Education of Students and Our Customers

We firmly believe that sharing our knowhow, and experience is very important to ensure successful and effective implantology methods and procedures. we provide hands on courses around the world, where we train graduate students and our customers to be as professional and educated\updated in the latest dental implantology procedures and innovative methods developed by our team of experts. We Adopt and Implement Digital Means and Technologies into The Dental World We embrace the technological changes involved in the digitalization of the dental world in order to support the present and future needs of our customers. Consequently, our digital CAD/CAM line offers a wide rangeof restoration products for our two implant connections. Additionally, our Guided Surgery Tool Kit supports guided surgery methods and enables dentists to select the software to use when planning surgeries, making their work easier and simpler.

We are environment friendly

We are very committed to environment keeping and take maximum measures to avoid using harmful and non-green materials in production lines we also encourage recycling of packing materials and prefer suppliers respecting environment. We encourage contributing to global community, We at Bioline participate in local-global events of contributing to community in developing countries and support pro bono contributor's by giving free products to aid needing communities.

R&d departments

Engineering cad-cam

Our company specializes in designing computer-aided dental elements. The planning department engineers have extensive experience in many areas.

Computerized mechanical design.

- Strength of materials.
- Simulation of assemblies and subassemblies.
- Metallurgy.

Human engineering and product planning to use a simple, stable, and comfortable to the doctor and the patient. Dental product design requires a comprehensive scientific thinking concerning mechanical strength of the product and the ability to ensure the survival during transplantation as well as the creation of a precise & strong connection and resistant to vibration and continuous wear. We are working on the continuous development of new products and improving existing products to ensure total quality so as to assure long-term success in improving the patient's quality of life.

State of the art accurate machining

Our company has a modern factory with high quality computerized machines, very secure, continuous, and repeatable accuracy. Manufacturing operations based on a computerized manufacturing system and program management and monitoring of the different jobs and very strict production process, controlled by skilled workers with knowledge and experience.

Quality inspection

Our company manufactures medical-dental products during application and adaptation to most stringent workmanship standards in the industry. The company's products are tested in the production process and in the transition between the various production stages by qualified inspectors in accordance with the specifications defined in production procedures. Our company is certified to ISO 13485-2012 and carries the CE-1023 icon. Our company is working constantly to improve the quality policy and ensuring maximum quality level of the company's products, all to ensure excellent products and customer satisfaction.

Clean room packing

Packaging process is done in a clean room, ensuring purification level specified by the medical device for dental implants. Packaging room environment ensures hygiene product during the sterile packaging and ensuring long shelf life without any possible contamination of the product. This room clean environment is suitable for international labor standards defined in the CE standard.



At Bioline, quality management is an integral part of the way we do business. Quality is one of our company's core values. We believe that in order to provide quality products and service to our customers, we must think-act-implement quality at all times. Our team is comprised of highly professional and qualified personnel who all contribute, each with their own expertise, to cooperatively create the high standards we uphold. Integrity and transparency are an integral part of our organizational culture, enabling us to offer products from the highest of standards in the dental market. Bioline dental Implant in-house manufacturing facilities allow it to efficiently produce commercial quantities of implants and accurate prosthetic parts to meet US, EU, Rus, India Mexico and many other countries requirements. We have a state of the art production facility that is dedicated to providing high quality products. Each and every one of our products designed and manufactured to meet the highest level of standards worldwide and we are regulatory approved by a number of leading notify bodies.

International Regulations-Certifications

As an advanced medical certified manufacturing facility, Bioline is fully compliance with all European directive for Medical Devise EN 93/42/EEC and ISO 13485:2016. Our factory is inspected on a regular basis by European Notified Body to ensure we are comply the requirements of the Standards.

Bioline is authorized to operate within USA and the EEA and all regulating bodies have given their approval to Bioline. We have invested a significant amount of resources to ensure that we maintain high quality in our manufacturing processes.

Bioline had gain and work according to the following standards:

CE marking - Certificate Number 14 0611 QS/NB/a

EN ISO 13485:2016 – Certificate Number I14081.

FDA listing of all prosthetic part and surgical tools


RUS – Russian regulation and marketing certificate

Free sale & Declarations of Conformity for medical device products are available

Lifetime Warranty Policy

Bioline has complete confidence in the quality and long-term success of the products we manufacture and market, therefor we offer a lifetime guarantee to all of our products.

provided they are used properly and in accordance with our instructions for use and maintenance.



Labeling & packaging



State of the art packaging, easy to use & to handle during the time of surgery.
clear warranty lable's contain the folowing simbles:

- BIO-GCM4510** ———• Product name + sizes
- | | |
|-----|--------------|
| REF | ILS1137 |
| LOT | BIO170290-01 |

 } ———• REF + LOT codes
- 2022-09 ———• expiration date
- | | |
|---------|---|
| STERILE | R |
|---------|---|

 ———• this product is sterilized by gamma radiation
- Do not Double sterilize
- for single use only



What is osseointegration?

Osseointegration, defined as a direct structural and functional connection between ordered, living bone and the surface of a load-carrying implant, is critical for implant stability, and is considered a prerequisite for implant loading and long-term clinical success of endosseous dental implants. The implant-tissue interface is an extremely dynamic region of interaction. This complex interaction involves not only biomaterial and biocompatibility issues but also alteration of mechanical environment. The processes of osseointegration involve an initial interlocking between alveolar bone and the implant body, and later, biological fixation through continuous bone apposition and remodeling toward the implant. The process itself is quite complex and there are many factors that influence the formation and maintenance of bone at the implant surface.

The usage & advantage of calcium phosphate:

Calcium phosphate (CaP) biomaterials currently in use for bone repair, substitution, augmentation, and regeneration include hydroxyapatite of synthetic or biologic origin, beta-tricalcium phosphate (β -TCP) and biphasic calcium phosphate. They are available as granules, porous blocks, CaP/polymer composites, cements, and as coatings on orthopedic and dental implants. Experimental CaP biomaterials include CO₃- and F-substituted apatites, Mg- and Zn-substituted β TCP, and CaP glasses.

Biofix : Biological Surface Treatment concept:

The use of a granular, multi-phase calcium phosphate abrasive such as **Biofix** Apatitic Abrasive has been recognized and Accepted for many years in the dental, orthopedic, and spinal markets.

Biofix Apatitic Abrasive provides a textured surface to either increase the surface area that comes in contact with the bone or as a surface preparation prior to the application of a coating such as hydroxyapatite or porous titanium.


Conventional abrasive blasting techniques that utilize alumina or silicon carbide always leave residual abrasive as a contaminant embedded in the implant surface.

Biofix Apatitic Abrasive will provide a residue-free surface after passivation.



The reference literature below links to research studies which discuss sandblasted and acid-etched (SLA) surface-treated implants—specifically with the use of Bionline's **Biofix** Apatitic Abrasive.

Our guiding line in **Biofix** formula is to avoid possible metal contamination caused by remains of abrasive materials sprayed on the titanium implant during surface treatment creation.

Biofix treatment provides the cleanest and most effective result to be monitored for many years to come.




Rotational implants group, 4-week healing period (H&E: A, $\times 40$; B, $\times 100$).




References:

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3602536/>
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3596629/>
http://www.jaypeejournals.com/eJournals/ShowText.aspx?ID=7317&Type=FREE&TYP=TOP&IN=_eJournals/images/JPLOGO.gif&IID=555&isPDF=YES


ESTHETIC ZONE RESTORATION USING SPIRAL DENTAL IMPLANTS By Prof. Dr Divye Malhotra



1. Pre-Op



2. Flap Opening



3. Osteotomy site preparation




4. Osteotomy site preparation



5. Implants placement's




6. Bone grafting + membrane



7. Bone grafting + membrane



8. suturing



9. X ray removed after 3 months – great integration



10. Abutments placement's



11. Final results

**CLINICAL & BIOLOGICAL INITIAL + LONG TERM STABILITY
,FAST HEALING, AND HIGH RATE OF OSSEOINTEGRATION
LEAD BY Biofix[®] FORMULA**

Drilling protocol

The recommendations presented in this drilling Protocol, based on long experience of specialists in the world. However it is important to examine each case on its merits based on the physician's ability to perform the procedure, and the state of the patient's jaw implant transpires with the beginning and there may be changes according to the situation.


In principle, we recommend using a drill that its diameter smaller than the diameter of the implant, in the case of soft bone condition. Also, consider using a similar drill diameter implant intended for implantation in the case of hard bone in order to avoid trauma to the bone and ensure rapid bone building around the implant.

IMPORTANT NOTICE :

It is important to ensure all instrumentation, surgical hand pieces, and equipment has been sterilized to prevent the possible contamination of the system, and the patient.

- recommended implant incretion torque 35-65 NCM
- recommended drilling speed 800-1200 RPM
- recommended prosthetics screw tightening torque 25-30 NCM
- recommended healing cap tightening torque 15 NCM
- all procedures recommended by BioLine Implants systems are not replacing the surgeon experience and judgment

Color Code





BioGuide – Guided surgery services

Currently, there have been major technological developments in diagnostic imaging methods and turned more accessible. The introduction of computed tomography images, including cone-beam computed tomography (CBCT), improved the outcome of implant-based treatments since it allows preoperative diagnosis with excellent quality and less exposure of the patients to radiation. When CBCT was associated with computer-aided design and manufacturing (CAD/CAM) provided surgical planning in both virtual and 3D environments, providing the practitioner with a realistic view of the patient's bony anatomy, thus permitting a virtual execution of the surgery in an ideal and precise prosthetically driven manner. This approach has been introduced to transfer virtual planning to the clinical procedure, allowing for less invasive surgeries, adequate implant placement, reduction of postoperative discomfort, and fabrication of prosthetic structures before surgical procedures.

The CAD/CAM has generated a growing trend in recent studies. There are different techniques involving guided surgery, evaluating accuracy of positioning of dental implants in comparison with virtual planning, and clinical results in the patients. The use of computer-guided surgery had been restricted to the advantage's surgical aspects of implant treatment. Prosthetic treatment still has to be carried out following conventional protocols. However, the link to transfer prosthetic information to the patient is of great importance, and exact reference points are required to position the implants in such a way that prefabricated prosthetics have a precise fit.

BioLine company established the BioGuide Service – computer guided surgery system in order to simplify implantology procedure to the beginners, and open the gate to the ones that was having a doubts. More than that BioGuide Service ensure high precision non-risk operation even for the experienced doctors, in a complicated case as implant placement in the esthetic zone / close to the mandibular nerve/ maxillary sinus / in a case of using extra-long implants as zygomatic / pterygoid as a graft less solution.

Contact us / your local dealer and plan your next safe & successful computer guided surgery!



Bio materials


Synthetic bioabsorbable polymers such as polylactic acid and polyglycolic acid used in various proportions, have been for many years, the most reliable components of numerous medical devices in the fields of dental, maxillofacial and orthopedic surgery, in the forms of synthetic; bone plates pins and screws, surgical sutures and shaped blades for maxillofacial applications. All these products have a high specific density principally because they must perform a mechanical action. Under similar implant conditions the absorption time is strictly related to the degree of density and can vary between 5-7 weeks to 2-3 years. Many variables are involved in biological degradation. These relate to the implant site, the age of the recipient, the tolerance of the tissue, as well as the physical aspects of the implant material, such as physical structure, chemical composition, molecular weight and form of exposed surface, and many others. The final by-products resulting from the metabolism of these polymers are carbon dioxide and water. Unlike products used in orthopedics and maxillofacial surgery which have an elevated mechanical resistance, SINTbone, based on a copolymer of polylactic and polyglycolic acid, is low in density, and is simply used as a space maintainer; and because of its low degree of density, it can guarantee absorption times between 4 to 8 months, in relation to the size of the implant, the reactivity of the recipient and the level of blood circulation at the implant site. SINTbone is used primarily as an absorbable biocompatible space maintainer. It is permeated by the blood and stabilizes the blood clot which is then progressively replaced by the osteoprogenitor cells to form new bone tissue. The spongy structure of SINTbone does not offer any resistance to the osteons and is physiologically absorbed which favors the conversion into bone tissue. Histological tests and SEM (Scansion Electron Microscopy) analysis have confirmed that the tissue does not suffer any unwanted short, medium or long-term reactions. All bio-compatibility tests carried out first on animals and then on humans have demonstrated that SINTbone is perfectly bio-compatible. SINTbone is available in three types: sintered BLOCK, sintered GRANULAR and SLURRY gel. The three different types make it easier to use SINTbone in relation to the type of surgery concerned and the site of the implant, because each type can be used either individually or in combination with the others.

Sintbone Granular Sintered



Unit weight : 450mg
Unit Volume : 1,6-1,8 ml ca.
Composition : Copolymer of
Polylactic and Polyglycolic Acid, Dextran

Sintbone Slurry Gel




Unit weight : 700mg
Unit Volume : 0,65 ml
Composition : Copolymer and homopolymers
of lactic and glycolic acid and, PEG, Dextran,
Distilled water Distilled water

Sintbone Block Sintered



Unit weight : 310mg
Unit Volume : 0,72 ml ca.
Composition : Copolymer of
Polylactic and Polyglycolic Acid, Dextran



INTERNAL HEX CONNECTION




Internal hex advantage's

- Easy simple to use
- One platform fits to all implant diameters
- Proven and well known
- High mechanical strength capacity
- High accuracy and easy future maintenance


SIMPLE & ACCURATE INTERNAL HEX COLLECTION




SDI




GRIP



MAI



CDI



PTR



ZYGO

Our spiral implant is a unique product based on traditional and proven design of tapered implant threaded balanced and dynamic self-drilling capability. Outstanding tapered implant symmetrical structure that contributes to balance efforts and equal distribution of possible loads during implantation and for the duration of use.


Anodized coating technology unique anti-bacterial process.

The implants undergo anodized coating process which contributes to creating anti-bacterial shell on titanium raw. The coating is made from natural ingredients Food Standard \ Medicine and unattended polluting industrial materials. Anodized great extent contributes to the mechanical strength and maintaining internal connection implant restoration parts and prevents friction between bare metal, this procedure ensures long-term success of rehabilitation on implants.

SDI Implant advantage's

- user friendly
- easy & fast incresion lead by the active design improved bone cutting and condensing ability's
- higher BIC (Bone & implant contact) lead for excellent primary stability
- can be used for immediate placement & loading reduce Bone resorption
- **Biofix** Biological surface for fast healing and long term osseointegration

RESTORATION OPTIONS:





Ordering information:

| | Dia ▶ | 2.9 | 3.5 | 3.75 | 4.2 | 5.0 | 6.0 |
|---------------|-------|------------|-----------|-----------|-----------|-----------|-----------|
| Length | 6 | | BIO-F3506 | BIO-F3706 | BIO-F4206 | BIO-F5006 | BIO-F6006 |
| | 8 | *BIO-F2908 | BIO-F3508 | BIO-F3708 | BIO-F4208 | BIO-F5008 | BIO-F6008 |
| | 10 | *BIO-F2910 | BIO-F3510 | BIO-F3710 | BIO-F4210 | BIO-F5010 | BIO-F6010 |
| | 11.5 | *BIO-F2911 | BIO-F3511 | BIO-F3711 | BIO-F4211 | BIO-F5011 | BIO-F6011 |
| | 13 | *BIO-F2913 | BIO-F3513 | BIO-F3713 | BIO-F4213 | BIO-F5013 | BIO-F6013 |
| | 16 | *BIO-F2916 | BIO-F3516 | BIO-F3716 | BIO-F4216 | BIO-F5016 | |

2.1 mm hex slim version

- The apical diameter is thinner than the Crestal diameter in 0.5 mm, included all implants sizes.


Spiral Implant



- Anti-bacterial anodized coating
- Micro rings -greater surface area
- Symmetrical neck – better load distribution



- Tapered body
- High condensing threads
- Sharp threads -self osteotomy creation
- High Primary Stability



- Sharp & deep thread for easy incision
- Narrow core - self tapping
- Prevent damage to the anatomical structures

- The Grip implant system has been designed especially to give better stability during immediate loading.
- Grip sharp & deep threads provide extremely strong retention for primary implant stability.
- The threads aggressive shape creates primary stability and maximum bone condensing in the sponge bone.
- Polished neck with groves prevents bone loss and increase healing and bonding with wider implant surface.
- The Grip implant designed to be used in extraction sites for immediate implantation.

Grip implant advantage's

- user friendly
- minimal drilling protocol required
- implant can be used as a self-osteotome
- design specially for under drilling sockets
- Gold anti-bacterial neck, prevent peri implantitis and crestal bone loss
- bone condensing threads-maximum primary implant stability
- **Biofix** Biological surface for fast healing and long term osseointegration


RESTORATION OPTIONS:



Ordering information:

| | Dia ▶ | 3.75 | 4.2 | 5.0 |
|--------|-------|--------------|--------------|--------------|
| Length | 8 | BIO-GRIP3708 | BIO-GRIP4208 | BIO-GRIP5008 |
| | 10 | BIO-GRIP3710 | BIO-GRIP4210 | BIO-GRIP5010 |
| | 11.5 | BIO-GRIP3711 | BIO-GRIP4211 | BIO-GRIP5011 |
| | 13 | BIO-GRIP3713 | BIO-GRIP4213 | BIO-GRIP5013 |
| | 16 | BIO-GRIP3716 | BIO-GRIP4216 | BIO-GRIP5016 |

Implant



Ø 4.0

3.0 mm

Ø 1.5

- Anti-bacterial anodized coating- bone loss prevention
- Micro rings -greater surface area
- Symmetrical neck – better load distribution
- Build in platform switching

- Extra tapered body
- High condensing threads
- Sharp threads -self osteotomy creation
- High Primary Stability

- Extra sharp & thin thread for easy incision
- Narrow core - self tapping
- Prevent damage to the anatomical structures

Our grip is a unique implant which is extra tapered to create extra self-osteotomy, bone condensing & preservation, by extra sharp and deep thread design+ thin and sharp apex (apex diameter is 1.5mm in all the implant diameters)

This unique featured design required different drilling sequence as per the Drill Less Approach

Please follow the sequence below:

IMPORTANT NOTICE :

It is important to ensure all instrumentation, surgical hand pieces, and equipment has been sterilized to prevent the possible contamination of the system, and the patient.


Soft bone (D3,D4): all sizes of grip implant to be drilled 2 mm less then the implant length

Hard bone (D2,D1): Ø2mm diameter (pilot drill) to be drilled full according to the implant length next drills to be drilled 2mm less then the implant length


final drill to be countersink (for 3.75&4.2 drill Ø3.75 - Ø4.2, for Ø5.0 implant drill Ø5-6)

- recommended implant incretion torque 35-65 NCM
- recommended drilling speed 800-1200 RPM
- recommended prosthetics screw tightening torque 25-30 NCM
- recommended healing cap tightening torque 15 NCM
- all procedures recommended by BioLine Implants systems are not replacing the surgeon experience and judgment

Color Code



drilling protocol



Active implant

is a unique dental implant developed after thorough research related to the ongoing quality of the connection between the bone and the dental implant. Our development team sought a solution that would extend as far as possible the life of the connection between the bone and the implant everything to ensure long-term quality restoration or almost infinite. In the development process is taken into consideration that even if we use the finest dental implant in the world, it is agreed that this object is a passive metal implant which has no lasting benefit or active contribution of osteointegration. from this point of view, we developed Active implant based on our spiral implant line, but hollowed and with few points of holes in the body, all for the transport or injection of fluid-gel that contain regeneration ingredients.

Inventive Thinking

As presented in the introduction about the Active family implant. The goal is to use this unique dental implant as a tool for transporting liquids and \ or special solutions designed to refresh the soft bone in the jaw. Special holes at the bottom of the implant body and graft are used for transporting liquids from the implant after successful implantation process and indefinitely when the implant is functioning as a basis for dental rehabilitation on the one hand, and as a tool for conveying the same fluid causes the processing of soft bone regeneration. This innovative Active implant is unique and exclusive to our company and reflects technological innovation and a major breakthrough in relation to prolong the life of dental implants and extending the medical care success and life quality of the patient. The proposed maintenance using Active dental implant \ device


RESTORATION OPTIONS:



Ordering information:

| | Dia ▶ | 3.3 | 3.75 | 4.2 | 5.0 | 6.0 |
|---------------|-------|------------|------------|------------|------------|------------|
| Length | 6 | | BIO-IR3706 | BIO-IR4206 | BIO-IR5006 | BIO-IR6006 |
| | 8 | BIO-IR3308 | BIO-IR3708 | BIO-IR4208 | BIO-IR5008 | BIO-IR6008 |
| | 10 | BIO-IR3310 | BIO-IR3710 | BIO-IR4210 | BIO-IR5010 | BIO-IR6010 |
| | 11.5 | BIO-IR3311 | BIO-IR3711 | BIO-IR4211 | BIO-IR5011 | BIO-IR6011 |
| | 13 | BIO-IR3313 | BIO-IR3713 | BIO-IR4213 | BIO-IR5013 | BIO-IR6013 |
| | 16 | BIO-IR3316 | BIO-IR3716 | BIO-IR4216 | BIO-IR5016 | |

Maintained Active Implant



- Anti-bacterial anodized coating
- Micro rings -greater surface area
- Symmetrical neck – better load distribution

- Innovative hollow platform ensure bone preservation & internal osseointegration with the years to come

- Sharp & deep thread for easy incision
- Narrow core - self tapping
- Prevent damage to the anatomical structures

The company's implant Cylinder built in order to achieve stability and serve as an anchor for a stable reconstruction. Implant corrugation delicate and multi-donor bone compression insertion process, the cylinder implant is recommended for use in cases of particularly hard bone, but functions well in other cases. cylinder Implant is a classic and reliable product for many years. Easy insertion of the implant in accordance with the recommended drilling Protocol and the geometric structure ensures long term stability.

Anodized coating technology unique anti-bacterial process.

The implants undergo anodized coating process which contributes to creating anti-bacterial shell on titanium raw. The coating is made from natural ingredients Food Standard \ Medicine and unattended polluting ontributes to the mechanical strength and maintaining internal connection implant restoration parts and prevents friction between bare metal, this procedure ensures long-term success of rehabilitation on implants.

CDI implant advantag's

- user friendly
- suit wide range of cases & bone types
- engineered for mild compression
- offer excellent stability in hard bone case's
- **Biofix** biological surface for fast healing and long term osseointegration

RESTORATION OPTIONS



Over Denture
straight/angulated till 30°


Screw Retained
straight/angulated till 50°

Cement Retained
straight/angulated till 45°


Ordering information:

| | Dia ▶ | 3.3 | 3.75 | 4.2 | 5.0 | 6.0 |
|--------|-------|-----------|-----------|-----------|-----------|-----------|
| Length | 6 | | BIO-G3706 | BIO-G4206 | BIO-G5006 | BIO-G6006 |
| | 8 | BIO-G3308 | BIO-G3708 | BIO-G4208 | BIO-G5008 | BIO-G6008 |
| | 10 | BIO-G3310 | BIO-G3710 | BIO-G4210 | BIO-G5010 | BIO-G6010 |
| | 11.5 | BIO-G3311 | BIO-G3711 | BIO-G4211 | BIO-G5011 | BIO-G6011 |
| | 13 | BIO-G3313 | BIO-G3713 | BIO-G4213 | BIO-G5013 | BIO-G6013 |
| | 16 | BIO-G3316 | BIO-G3716 | BIO-G4216 | BIO-G5016 | |

Cylindrical Implant



- Anti-bacterial anodized coating
- Micro rings -greater surface area
- Symmetrical neck – better load distribution



- Straight body
- condensing threads
- High Primary Stability for d1/d2 type of bone



- Deep condensing grooves
- Stable core apex
Prevent damage to the anatomical structures

This implant is a single unit containing a body comprising an implant and a straight structure.

The implant is made of 3 mm diameter 3.3 mm and 3.75 mm.

This implant is designed specifically for use in a narrow ridge and narrow place between teeth.

The implant is for immediate loading because of its uniform structure, this implant coated with anti-bacterial anodized gold colored.

The surface of the implant treated the same way in which the treatment is carried out at all our implants series.


OPI Implant advantage

- build in platform switching for soft tissue growth
- user friendly especially for narrow ridges and tight spaces
- can be used for immediate loading
- smooth anti-bacterial collar prevents peri implantitis
- self-tapping for easy incision
- **Biofix** Biological surface for fast healing and long term osseointegration

Ordering information:

| | Dia ▶ | 3.0 | 3.3 | 3.75 |
|--------|-------|-----------|-----------|-----------|
| Length | 10 | BIO-P3010 | BIO-P3310 | BIO-P3710 |
| | 11.5 | BIO-P3011 | BIO-P3311 | BIO-P3711 |
| | 13 | BIO-P3013 | BIO-P3313 | BIO-P3713 |
| | 16 | BIO-P3016 | BIO-P3316 | BIO-P3716 |


One Piece Implant



- Anti-bacterial anodized coating
- Single piece – build in abutment
- Symmetrical neck – better load distribution


- Tapered body
- High condensing threads
- Sharp threads -self osteotomy creation
- High Primary Stability

- Sharp thread for easy incision
- Narrow core - self tapping
- Prevent damage to the anatomical structures




The TILT Concept

**Innovative GraftLess Solution For
Atrophic Jaws Using Immediate
Loading Technique.**



Zygomatic implants are an evidence-based surgical and prosthetic solution for both two-stage and immediate loading protocols. Today, zygomatic implants are usually placed using an immediate loading protocol. The main indication for zygomatic implants is the severely resorbed edentulous maxilla, but they can also be used in partially edentulous situations. Indications for zygomatic implant insertion include: alternative for sinus augmentation, failed sinus augmentation, rehabilitation after tumor resection or trauma, failure of conventional implants, failure of previous bone grafts. The placement of zygomatic implants requires adequate training and surgical experience.



Zygomatic Implant advantage

- Graftless solution – avoid sinus augmentations-time & money saver
- Smooth anti-bacterial long neck – prevent inflammation, even in intra sinus cases
- Shorten the time of permanent teeth delivery
- high primary stability as they supported by the zygomatic bone
- high patient satisfaction

RESTORATION OPTIONS:



Screw Retained
straight/angulated till 50°


Cement Retained
straight/angulated till 45°

Ordering information:

| | Dia ▶ | 3.75 | 4.2 |
|--------|-------|--------------|--------------|
| Length | 30 | BIO-ARCH3730 | BIO-ARCH4230 |
| | 35 | BIO-ARCH3735 | BIO-ARCH4235 |
| | 40 | BIO-ARCH3740 | BIO-ARCH4240 |
| | 45 | BIO-ARCH3745 | BIO-ARCH4245 |
| | 50 | BIO-ARCH3750 | BIO-ARCH4250 |
| | 55 | BIO-ARCH3755 | BIO-ARCH4255 |
| | 60 | BIO-ARCH3760 | BIO-ARCH4260 |

Important notice: 3.75 mm diameter of zygomatic implant to be used in case of quad maxillary zygoma and should be placed beside 4.2 mm to split the load distribution and avoid implants fractures

Zygomatic Implant



- Anti-bacterial anodized coating
- Symmetrical neck – better load distribution

- Smooth long neck – prevent bacteria & inflammation

- Sharp & deep thread for easy incision
- Narrow core - self tapping
- 13mm threaded body – perfect fits to the zygoma bone

Pterygoid implant insertion is an alternative to avoid sinus-lifting or other grafting procedures to treat the posterior maxilla

Pterygoid implant are especially used in partial edentulism in order to avoid distal cantilevers

The placement of a pterygoid implant requires surgical experience

Pterygoid implants have high success rates, lesser bone loss levels to those of conventional implants, minimal complications and a good patient acceptance

Pterygoid implant advantage's

- Graftless solution – avoid sinus augmentations-time & money saver
- Shorten the time of permanent teeth delivery
- Gold anti-bacterial neck, prevent peri implantitis and crestal bone loss
- bone condensing threads-maximum primary implant stability
- Sharp thin apex, easy biocritical anchorage
- **Biofix** Biological surface for fast healing and long term osseointegration

RESTORATION OPTIONS:



Ordering information:

| | Dia ▶ | 3.0 | 3.5 | 3.75 | 4.2 |
|---------------|-------|----------------------------------|-----------|-----------|-----------|
| Length | 16 | BIO-I3016 | BIO-I3516 | BIO-I3716 | BIO-I4216 |
| | 18 | BIO-I3018 | BIO-I3518 | BIO-I3718 | BIO-I4218 |
| | 20 | BIO-I3020 | BIO-I3520 | BIO-I3720 | BIO-I4220 |
| | 22 | BIO-I3022 | BIO-I3522 | BIO-I3722 | BIO-I4222 |
| | 25 | BIO-I3025 | BIO-I3525 | BIO-I3725 | BIO-I4225 |
| | | * Slim Version 2.1 Mm Hexagon | | | |

Pterygoid Implant



- gold anodized anti bacterial coating
- machined micro thread design
- crestal bone preservation
- high initial stability



- high cutting
- high bone condensing aggressive
- square threaded design



- sharp apex for bicortical anchorage
- gold anodized anti bacterial coating
- easy insertion
- high primary engagement

The concept of tilted implants brings to the table a unique state of the art prosthetics components line in order to improve the quality of restoration which was the only limit of tilted implants many years back. The use of our innovative abutment's makes restoration easier, with the perfect achievement of precision and correct the angulation of the implant after placement without a limit. Our all premium tilted implants prosthetics supported the platform switching concept with a thin platform level available in different heights to ensure perfect fitting in sub - crestal implant placement that prevent crestal bone loss and increase the success rate of the implant even in compromise situations

A Revolution In Dental Implantology





one smart abutment

Straight Abutment

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-IS-05 | BIO-IS-15 | BIO-IS-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |



15° Angulated Abutment

| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-IA15-05 | BIO-IA15-15 | BIO-IA15-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |



25° Angulated Abutment

| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-IA25-05 | BIO-IA25-15 | BIO-IA25-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |



35° Angulated Abutment

| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-IA35-05 | BIO-IA35-15 | BIO-IA35-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |



45° Angulated Abutment

| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-IA45-05 | BIO-IA45-15 | BIO-IA45-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |




Angulated Multi-Unit Abutment

| | | | | | |
|---------|--------------|--------------|--------------|--------------|--------------|
| Cat No. | BIO-MS-09-10 | BIO-MS-18-10 | BIO-MS-30-10 | BIO-MS-40-10 | BIO-MS-50-10 |
|---------|--------------|--------------|--------------|--------------|--------------|

* All angles are available height 1mm

**IMMEDIATE IMPLANTATION LOADING WITHOUT CANTILEVER
USING SCREW RETAINED PROSTHETICS
MAXILLARY ARCH REHABILITATION
By Dr P. Venkat Nag. Hyderabad India**



1. Maxillary arch failing dentition lower 36,46,47 loeer implants placed 2 years back




2. Pre operation




3. Single drill concept osseointegration mult cortical anchorage



4. Biloline pterygoid (implant) sharp apex for easy cortical engagement



5. Nasal cortical bone engagement



6. Pterygoid cortical bone engagement



7. Reverse torque test >45 NC




8. High torque >60 NC bicortical engagement




9. Pterygoid implant before placement



10. Flapless pterygoid engaging the pterygoid pillar




11. All on 6 pterygoid implants




12. Impressions done on the day of surgery



13. Final results



14. Immediate loading + no cantilever + screw retained prosthesis in just 2 days with all tilt protocol





nerve bypass technique


INTERNAL HEX

prosthetics
solutions

healing abutments



angulated abutment



multi unit

Healing Caps

Diameter 3.3 mm Slim *for 2.9 mm implant

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-H3302 | BIO-H3303 | BIO-H3304 | BIO-H3305 | BIO-H3306 | BIO-H3307 |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |



Diameter 3.8 mm Standard

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-H3802 | BIO-H3803 | BIO-H3804 | BIO-H3805 | BIO-H3806 | BIO-H3807 |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |



Diameter 4.6 mm Standard

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-H4602 | BIO-H4603 | BIO-H4604 | BIO-H4605 | BIO-H4606 | BIO-H4607 |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |




Diameter 5.5 mm Wide

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-H5502 | BIO-H5503 | BIO-H5504 | BIO-H5505 | BIO-H5506 | BIO-H5507 |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |



Diameter 6 mm Wide

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-H6302 | BIO-H6303 | BIO-H6304 | BIO-H6305 | BIO-H6306 | BIO-H6307 |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |



Healing
Caps

Straight Abutments



Straight Abutment

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-A5908 | BIO-A5909 | BIO-A5911 | BIO-A5912 | BIO-A5601 | BIO-A5602 |
| Height | 8.5 mm | 9.5 mm | 10.5 mm | 11 mm | 8.5 mm | 10.5 mm |



Straight Shoulder Abutment

| | | | | | | |
|----------|-----------|-----------|-----------|-----------|-----------|----------------------|
| Cat No. | BIO-A4801 | BIO-A4802 | BIO-A4803 | BIO-A4804 | BIO-A5001 | BIO-A5002 Non-Hex |
| Height | 8.5 mm | 9.5 mm | 10.5 mm | 11.5 mm | 9.5 mm | 9.5 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm | 1.70 mm | 1.70 mm |



Straight Shoulder Abutment

| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-A5300 | BIO-A5301 | BIO-A5302 |
| Height | 8.5 mm | 8.5 mm | 8.5 mm |
| Shoulder | 0.50 mm | 1.50 mm | 2.50 mm |



Straight smooth Abutment

| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-A6006 | BIO-A6008 | BIO-A6011 |
| Height | 6.00 mm | 8.00 mm | 11.00 mm |
| Shoulder | N/A | N/A | N/A |



Straight Anatomic Abutment

| | | | | |
|----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-A4601 | BIO-A4602 | BIO-A4603 | BIO-A4604 |
| Height | 9.00 mm | 10.00 mm | 11.00 mm | 12.00 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm |



Slim straight Abutment *for 2.9 mm implant

| | | |
|---------|-------------|-------------|
| Cat No. | BIO-A5909-S | BIO-A5911-S |
| Height | 9.0 mm | 11.0 mm |



Slim straight Anatomic Abutment *for 2.9 mm implant


| | | | | |
|----------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-A3801-S | BIO-A3802-S | BIO-A3803-S | BIO-A3804-S |
| Height | 9.00 mm | 10.00 mm | 11.00 mm | 12.00 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm |



Screw for Abutment

Included with all types of abutments

| | |
|---------|-----------|
| Cat No. | BIO-S8324 |
|---------|-----------|



Straight
Abutments

Angulated Abutments

15°

15° Standard Abutment

| | | | | |
|---------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-A3207 | BIO-A3208 | BIO-A3209 | BIO-A3211 |
| Height | 7 mm | 8.5 mm | 9 mm | 11 mm |



15° Anatomic Standard Abutment

| | | | | | |
|----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-A2400 | BIO-A2401 | BIO-A2402 | BIO-A2403 | BIO-H3806 |
| Height | 8.5 mm | 9.5 mm | 10.5 mm | 11.5 mm | 12.5 mm |
| Shoulder | 0.80 mm | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm |



15° Standard Slim Abutment *for 2.9 mm implant

| | | |
|---------|-------------|-------------|
| Cat No. | BIO-A3209-S | BIO-A3211-S |
| Height | 9 mm | 11 mm |



15° Anatomic Slim Abutment *for 2.9 mm implant

| | | | | |
|----------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-A2401-S | BIO-A2402-S | BIO-A2403-S | BIO-A2404-S |
| Height | 9.5 mm | 10.5 mm | 11.5 mm | 12.5 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm |



25°

25° Standard Abutment

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-A3407 | BIO-A3409 | BIO-A3411 |
| Height | 7 mm | 9 mm | 11 mm |



25° Standard Slim Abutment *for 2.9 mm implant

| | |
|---------|-------------|
| Cat No. | BIO-A3407-S |
| Height | 7 mm |



25° Anatomic Standard Abutment

| | | | | |
|----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-A2601 | BIO-A2602 | BIO-A2603 | BIO-A2604 |
| Height | 9.5 mm | 10.5 mm | 11.5 mm | 12.5 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm |



35°

35° Standard Abutment


| | |
|---------|-----------|
| Cat No. | BIO-A3509 |
| Height | 9 mm |



45°

45° Standard Abutment

| | |
|---------|-----------|
| Cat No. | BIO-A4509 |
| Height | 9 mm |



Angulated Abutments

Temporary Peek Abutments

Peek Straight Anatomic Abutment



| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-P4510 | BIO-P4520 | BIO-P4530 |
| Height | 9 mm | 10 mm | 11 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm |

Peek Anatomic Abutment 15°



| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-P2401 | BIO-P2402 | BIO-P2403 |
| Height | 9 mm | 10 mm | 11 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm |

Peek Anatomic Abutment 25°



| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-P2601 | BIO-P2602 | BIO-P2603 |
| Height | 9 mm | 10 mm | 11 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm |

Castable Abutments



Castable Narrow Abutment

| | | |
|---------|-----------|-----------|
| Cat No. | BIO-C1001 | BIO-C1002 |
| | NON-HEX | HEX |



Slim Castable Abutment *for 2.9 mm implant

| | | |
|---------|-------------|-------------|
| Cat No. | BIO-C1001-S | BIO-C1002-S |
| | NON-HEX | HEX |



Castable Abutment

| | | |
|---------|-----------|-----------|
| Cat No. | BIO-C2003 | BIO-C2004 |
| | NON-HEX | HEX |



Ucla Abutment


| | |
|---------|-----------|
| Cat No. | BIO-C2001 |
| | UCLA |



Screw for Abutment

Included with all types of abutments

| | |
|---------|-----------|
| Cat No. | BIO-S8324 |
|---------|-----------|



Temporary Peek Abutments

Impression components

Transfer For Open Tray

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-T4008 | BIO-T4201 | BIO-T4012 |
| Height | 4.75 mm | 4.30 mm | 4.75 mm |
| | 8.00 mm | 12.00 mm | 12.00 mm |



Slim Transfer For Open Tray *for 2.9 mm implant

| | |
|---------|-------------|
| Cat No. | BIO-T4012-S |
| Height | 4.75 mm |
| | 12.00 mm |



Transfer For Close Tray

| | | | |
|----------|-----------|-----------|-----------|
| Cat No. | BIO-T3507 | BIO-T3511 | BIO-T3601 |
| Diameter | 3.85 mm | 4.50 mm | 4.75 mm |
| Height | 12.00 mm | 12.00 mm | 8.00 mm |



Slim Transfer For Close Tray *for 2.9 mm implant

| | | |
|----------|-------------|-------------|
| Cat No. | BIO-T3507-S | BIO-T3601-S |
| Diameter | 3.85 mm | 4.75 mm |
| Height | 12.00 mm | 8.00 mm |



Click Transfer For Close Tray

| | | |
|----------|-----------|-----------|
| Cat No. | BIO-T3409 | BIO-T3413 |
| Diameter | 4.70 mm | 4.70 mm |
| Height | 9.00 mm | 13.00 mm |



Ball Transfer

| | | | | | |
|----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-T3801 | BIO-T3802 | BIO-T3803 | BIO-T3804 | BIO-T3805 |
| Diameter | 4.75 mm | 4.75 mm | 4.75 mm | 4.75 mm | 4.75 mm |
| Shoulder | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm |



Plastic Cap For Ball Transfer

| | | |
|----------|-----------|-----------|
| Cat No. | BIO-T4401 | BIO-T4402 |
| Diameter | 4.75 mm | 4.75 mm |
| Height | 7.00 mm | 10.00 mm |



Screw for Transfer

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-S1307 | BIO-S1610 | BIO-S2418 |
| Height | 13.00 mm | 16.00 mm | 24.00 mm |




Implant Analog

| | | |
|----------|-----------|-----------|
| Cat No. | BIO-T6404 | BIO-T6405 |
| Diameter | 4.00 mm | 5.00 mm |
| Height | 12.75 mm | 12.75 mm |



Slim Implant Analog *for 2.9 mm implant

| | |
|----------|-------------|
| Cat No. | BIO-T6404-S |
| Diameter | 4.00 mm |
| Height | 12.75 mm |



Impression Components

Ball attachment



Ball Attachment

| | | | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-T1200 | BIO-T1201 | BIO-T1202 | BIO-T1203 | BIO-T1204 | BIO-T1205 | BIO-T1206 | BIO-T1207 |
| Height | 0.5 mm | 1 mm | 2 mm | 3 mm | 4 mm | 5 mm | 6 mm | 7 mm |



Slim Ball Attachment *for 2.9 mm implant

| | | | | | | | | |
|---------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| Cat No. | BIO-T1200-SL | BIO-T1201-SL | BIO-T1202-SL | BIO-T1203-SL | BIO-T1204-SL | BIO-T1205-SL | BIO-T1206-SL | BIO-T1207-SL |
| Height | 0.5 mm | 1 mm | 2 mm | 3 mm | 4 mm | 5 mm | 6 mm | 7 mm |




Metal Cap For Ball Attachment

| | |
|---------|-----------|
| Cat No. | BIO-T3001 |
|---------|-----------|

Caps 2.5mm for Ball Attachment



| | | | |
|---------|------------------------|------------------|----------------------|
| Cat No. | BIO-T3002 | BIO-T3003 | BIO-T3004 |
| | Silicone Extra soft | Silicone Soft | Silicone Standart |



Ball Attachment

Locator attachment

Power Click Straight Abutment

| | | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-PCA00 | BIO-PCA01 | BIO-PCA02 | BIO-PCA03 | BIO-PCA04 | BIO-PCA05 | BIO-PCA06 |
| Height | 0.5 mm | 1 mm | 2 mm | 3 mm | 4 mm | 5 mm | 6 mm |



Metal Cap For Power Click Attachment

| | |
|---------|-----------|
| Cat No. | BIO-PCMC0 |
|---------|-----------|



Silicone Sealed Ring




| | |
|---------|-----------|
| Cat No. | BIO-PCSSR |
|---------|-----------|



Silicone Softness



| | | | | |
|---------|------------|-----------|-----------|-----------|
| Cat No. | BIO-CN088 | BIO-CN087 | BIO-CN086 | BIO-CN090 |
| | Extra soft | soft | standart | Hard |



Locator Attachment

Angulated attachment

Versatile Angulated Connection Base



| | | | |
|----------|-------------|-------------|-------------|
| Cat No. | BIO-VAC09-I | BIO-VAC18-I | BIO-VAC30-I |
| Diameter | 4.60 mm | 4.60 mm | 4.60 mm |
| Angle | 9° | 18° | 30° |

Ball Attachment For V.A.C




| | | | | |
|---------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-BVAC2-I | BIO-BVAC3-I | BIO-BVAC4-I | BIO-BVAC5-I |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm |

Power Click Attachment For V.A.C



| | | | | |
|---------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-PVAC2-I | BIO-PVAC3-I | BIO-PVAC4-I | BIO-PVAC5-I |
| Height | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm |



Angulated Attachment

Multi-units abutments

Multi-Unit Straight Abutment

| | | | | | | |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-M0101 | BIO-M0102 | BIO-M0103 | BIO-M0104 | BIO-M0105 | BIO-M0106 |
| Height | 1.00 mm | 2.00 mm | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm |



Slim Multi-Unit Straight Abutment *for 2.9 mm implant

| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-M0101-S | BIO-M0102-S | BIO-M0103-S |
| Height | 1.00 mm | 2.00 mm | 3.00 mm |



Multi-Unit Plastic castable Sleeve

| | |
|---------|-----------|
| Cat No. | BIO-MPL01 |
|---------|-----------|



Multi-Unit Titanium Sleeve

| | |
|---------|-----------|
| Cat No. | BIO-MTI01 |
|---------|-----------|



Multi-Unit Short Titanium Sleeve

| | |
|---------|-----------|
| Cat No. | BIO-MTI0S |
|---------|-----------|



Healing Cap For Multi-Unit

| | |
|---------|-----------|
| Cat No. | BIO-M0HC1 |
|---------|-----------|



Transfer For Multi-Unit

| | |
|---------|-----------|
| Cat No. | BIO-M0TR1 |
|---------|-----------|



Analog For Multi-Unit


| | |
|---------|-----------|
| Cat No. | BIO-M0AN1 |
|---------|-----------|



Screw for Multi-Unit

Included with all types of abutments

| | |
|---------|-----------|
| Cat No. | BIO-MOSC1 |
|---------|-----------|



Multi-Units
Abutments

Angulated multi-unit



Angulated Multi-Unit Abutment

Cat No. | BIO-MU09 | BIO-MU18 | BIO-MU30 | BIO-MU40 | BIO-MU50 |

* All angles are available heights 1-3



Slim Angulated Multi-Unit Abutment *for 2.9 mm implant

Cat No. | BIO-MU09-S | BIO-MU18-S | BIO-MU30-S | BIO-MU40-S | BIO-MU50-S |

* All angles are available heights 1-3



Multi-Unit Plastic Castable Sleeve

Cat No. | BIO-MPLD |



Multi-Unit Titanium Sleeve

Cat No. | BIO-MUTSDL |



Multi-Unit Short Titanium Sleeve

Cat No. | BIO-MUTSDH |



Healing Cap For Multi-Unit

Cat No. | BIO-MUHC1 |



Close Transfer For Multi-Unit

Cat No. | BIO-MUTRIC |



Open Transfer For Multi-Unit

Cat No. | BIO-MUTRIO |



Analog For Multi-Unit


Cat No. | BIO-MUANAD |




Screw for Angulated Multi-Unit

Included with all types of abutments

Cat No. | BIO-MSC |



Angulated
Multi-Unit



All the needs for a successful screw retained restoration in one user friendly premium set. Simplify the whole procedure for single / multiple units.


Multi-unit premium set contain:

- Straight / angulated multi-unit abutment
- Castable sleeve
- Titanium sleeve
- Healing cap
- Transfer copy
- Lab analog
- Prosthetic's screw



full arch impression coping
screw retained prosthetics



CONICAL CONNECTION IMPLANTS



Conical connection advantage's

- Easy simple to use
- Sealed connection reduce bacteria risk
- Prevent screw loosening
- High mechanical strength capacity
- High accuracy and easy future maintenance

NEW GENERATION OF CONICAL CONNECTION IMPLANT SERIES




CCR

Our series of conical connection implants based on our proved and familiar structure of spiral implant [SDI] which become a best seller proven product over more than 32 countries worldwide more than 12 years. Conical connection design brings the solution for ensuring long term stability of restoration parts connected to the implant. The conical geometry of implants connection area and prosthetic parts spread on wider basis and create cold welding between the parts with almost no load on the traditional titanium screw, this geometry also double the mechanical stance of the implant-abutment during the use by patient for many years to come. The connection point is 100% sealed and tightened.

CCR Implant advantage's

- user friendly
- easy & fast incision lead by the active design
- improved bone cutting and condensing ability's
- higher BIC (Bone & implant contact) lead for excellent primary stability
- can be used for immediate placement & loading
- reduce Bone resorption
- **Biofix** Biological surface for fast healing and long term osseointegration.

RESTORATION OPTIONS:



Ordering information:

| | Dia ▶ | 3.75 | 4.2 | 5.0 | 6.0 |
|---------------|-------|--------------|--------------|--------------|---------------|
| Length | 6 | *BIO-CCR3706 | *BIO-CCR4206 | *BIO-CCR5006 | *BIO-CCR6006 |
| | 8 | *BIO-CCR3708 | *BIO-CCR4208 | *BIO-CCR5008 | *BIO-CCR6008 |
| | 10 | *BIO-CCR3710 | *BIO-CCR4210 | *BIO-CCR5010 | *BIO-CCR6010 |
| | 11.5 | *BIO-CCR3711 | *BIO-CCR4211 | *BIO-CCR5011 | *BIO-CCR6011 |
| | 13 | *BIO-CCR3713 | *BIO-CCR4213 | *BIO-CCR5013 | *BIO-CCR6013 |
| | 16 | *BIO-CCR3716 | *BIO-CCR4216 | *BIO-CCR5016 | *BIO -CCR6016 |

Conical connection



- Anti-bacterial anodized coating
- Micro rings -greater surface area
- Symmetrical neck – better load distribution





- Tapered body
- High condensing threads
- Sharp threads -self osteotomy creation
- High Primary Stability




- Sharp & deep thread for easy incision
- Narrow core - self tapping
- Prevent damage to the anatomical structures

CONICAL CONNECTION
prosthetic's solutions

healing abutments



angulated abutment



multi unit

Conical Healing Cap

Diameter 4.0 mm Standard



| | | | | | |
|---------|-------------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-CHC4003 | BIO-CHC4004 | BIO-CHC4005 | BIO-CHC4006 | BIO-CHC4007 |
| Height | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |

Diameter 5.0 mm Wide



| | | | | | |
|---------|-------------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-CHC5003 | BIO-CHC5004 | BIO-CHC5005 | BIO-CHC5006 | BIO-CHC5007 |
| Height | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |

Diameter 6.0 mm Wide



| | | | | | |
|---------|-------------|-------------|-------------|-------------|-------------|
| Cat No. | BIO-CHC6003 | BIO-CHC6004 | BIO-CHC6005 | BIO-CHC6006 | BIO-CHC6007 |
| Height | 3.00 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |

Conical Straight Abutment



| | | | |
|---------|-------------|-------------|-------------|
| Cat No. | BIO-CCA5907 | BIO-CCA5909 | BIO-CCA5911 |
| Height | 7.00 mm | 9.00 mm | 11.00 mm |

Conical Anguletad Abutment 15°



| | | | | |
|---------|------------|------------|------------|------------|
| Cat No. | BIO-CA3207 | BIO-CA3208 | BIO-CA3209 | BIO-CA3211 |
| Height | 7.00 mm | 8.00 mm | 9.00 mm | 11.00 mm |

Conical Anguletad Abutment 25°



| | | | | |
|---------|------------|------------|------------|------------|
| Cat No. | BIO-CA3407 | BIO-CA3408 | BIO-CA3409 | BIO-CA3411 |
| Height | 7.00 mm | 8.00 mm | 9.00 mm | 11.00 mm |

C-C prosthetic's

Conical Ball Attachment

| | | | | | | | |
|---------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| Cat No. | BIO-CBA-1200-1-S | BIO-CBA-1200-2-S | BIO-CBA-1200-3-S | BIO-CBA-1200-4-S | BIO-CBA-1200-5-S | BIO-CBA-1200-6-S | BIO-CBA-1200-7-S |
| Height | 1.00 mm | 2.00 mm | 3.0 mm | 4.00 mm | 5.00 mm | 6.00 mm | 7.00 mm |



Conical Locator

| | | | | | | |
|---------|--------------|--------------|--------------|--------------|--------------|--------------|
| Cat No. | BIO-CCPCA-01 | BIO-CCPCA-02 | BIO-CCPCA-03 | BIO-CCPCA-04 | BIO-CCPCA-05 | BIO-CCPCA-06 |
| Height | 1.00 mm | 2.00 mm | 3.0 mm | 4.00 mm | 5.00 mm | 6.00 mm |



Conical straight Multi Unit

| | | | | | | |
|---------|------------|------------|------------|------------|------------|------------|
| Cat No. | BIO-CCM101 | BIO-CCM102 | BIO-CCM103 | BIO-CCM104 | BIO-CCM105 | BIO-CCM106 |
| Height | 1.00 mm | 2.00 mm | 3.0 mm | 4.00 mm | 5.00 mm | 6.00 mm |



Conical Angulated multi Unit

| | | | | | |
|---------|------------|------------|------------|------------|------------|
| Cat No. | BIO-CAMU09 | BIO-CAMU18 | BIO-CAMU30 | BIO-CAMU40 | BIO-CAMU50 |
| degrees | 9° | 18° | 30° | 40° | 50° |



Conical Transfer

| | | | |
|----------|------------|------------|------------|
| Cat No. | BIO-T3507C | BIO-T3601C | BIO-T4212C |
| diameter | 3.85 mm | 4.75 mm | 4.80 mm |
| height | 12.0mm | 8.0 mm | 14.8 mm |



Conical Lab Analog

| | |
|---------|-------------|
| Cat No. | BIO-T6404CC |
| height | 12.75 mm |



C-C prosthetic's



Conical Motor Mount

| | |
|---------|------------|
| Cat No. | BIO-X1015C |
| height | 28.0 mm |



Conical Implants Driver

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-X2307 | BIO-X2310 | BIO-X2315 |
| height | 7.0 mm | 10.0 mm | 15.0 mm |




Conical Abutment Remover


| | |
|---------|----------|
| Cat No. | BIO-HCRA |
| height | 7.0 mm |

**UNIVERSAL
INSTRUMENTS**
line for all
procedures

3 step drill



universal drill



implant driver

Drills



Countersink

| | | |
|----------|----------------|----------------|
| Cat No. | BIO-D1034 | BIO-D1056 |
| Diameter | ∅ 3.80-4.20 mm | ∅ 5.00-6.00 mm |



Marking

| | |
|---------|-----------|
| Cat No. | BIO-D3410 |
|---------|-----------|



Extension

| | |
|---------|-----------|
| Cat No. | BIO-D3412 |
|---------|-----------|



Standard 16 mm

| | | | | | | | |
|----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-D1220 | BIO-D1225 | BIO-D1228 | BIO-D1232 | BIO-D1236 | BIO-D1242 | BIO-D1252 |
| Diameter | ∅ 2.00 mm | ∅ 2.50 mm | ∅ 2.8 mm | ∅ 3.2 mm | ∅ 3.65 mm | ∅ 4.20 mm | ∅ 5.20 mm |



Trephine

| | | | | |
|----------|----------------|----------------|----------------|----------------|
| Cat No. | BIO-D2020 | BIO-D2030 | BIO-D2040 | BIO-D2050 |
| Diameter | ∅ 2.00-3.00 mm | ∅ 3.00-4.00 mm | ∅ 4.00-5.00 mm | ∅ 5.00-6.00 mm |



3 Level

| | | |
|---------|-----------|-----------|
| Cat No. | BIO-TDL11 | BIO-TDL18 |
| Length | 11.50 mm | 18.00 mm |

Ratchet

| Cat No. | BIO-X1020 |



Torque Ratchet

| Cat No. | BIO-X1021 |



Surgical Screw Driver Hex 6.35 mm

| Cat No. | BIO-X1023 |



Surgical Screw Driver Hex 2.42 mm

| Cat No. | BIO-X1022 |



Depth Probe

| Cat No. | BIO-X1025 |





Driver Hex 1.25 mm

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-X2107 | BIO-X1210 | BIO-X1215 |
| length | 7.0 mm | 10.0 mm | 15.0 mm |



Driver Hex 2.42 mm

| | | | | |
|---------|-----------|-----------|-----------|-----------|
| Cat No. | BIO-X2407 | BIO-X2410 | BIO-X2415 | BIO-X2423 |
| length | 7.0 mm | 10.0 mm | 15.0 mm | 23 mm |



Driver Hex 2.42 mm Front Ring

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-X2607 | BIO-X2610 | BIO-X2615 |
| length | 7.0 mm | 10.0 mm | 15.0 mm |



Slim Hex Driver 2.1 mm *for 2.9 mm implant

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-X2807 | BIO-X2810 | BIO-X2815 |
| length | 7.0 mm | 10.0 mm | 15.0 mm |



Driver Hex 1.25 mm long

| | |
|---------|-----------|
| Cat No. | BIO-X1006 |
| length | 15.0 mm |



Driver Hex 1.25 mm short

| | |
|---------|-----------|
| Cat No. | BIO-X1007 |
| length | 7.0 mm |



Starter Surgical Kit Contain:

- Rachet 6.35 hex
- Drill extender
- marking drill
- 2 pilot drill
- 2.8 drill
- 3.2 drill
- 3.65 drill
- 4.2 drill
- 5.2 drill
- 1.25mm hex prosthetic's driver
- 2.42 mm hex implant driver

- BIO-X1020
- BIO-D3412
- BIO-D3410
- BIO-D1220
- BIO D1228
- BIO D1232
- BIO-D1236
- BIO-D1242
- BIO-D1252
- BIO-X1210
- BIO-X2610



Safe kit



The Integral Stopper Drills Kit assures accurate depth control. The Kit is designed for each Implant Diameter and Length. Simplify drilling procedure!

Safe drills surgical kit contain:

| | | | |
|---------------|-------------|---------------|--------------|
| stopper drill | ∅ 2.0 L6 | stopper drill | ∅ 3.65 L6 |
| stopper drill | ∅ 2.0 L8 | stopper drill | ∅ 3.65 L8 |
| stopper drill | ∅ 2.0 L10 | stopper drill | ∅ 3.65 L10 |
| stopper drill | ∅ 2.0 L11.5 | stopper drill | ∅ 3.65 L11.5 |
| stopper drill | ∅ 2.0 L13 | stopper drill | ∅ 3.65 L13 |
| stopper drill | ∅ 2.8 L6 | stopper drill | ∅ 4.2 L6 |
| stopper drill | ∅ 2.8 L8 | stopper drill | ∅ 4.2 L8 |
| stopper drill | ∅ 2.8 L10 | stopper drill | ∅ 4.2 L10 |
| stopper drill | ∅ 2.8 L11.5 | stopper drill | ∅ 4.2 L11.5 |
| stopper drill | ∅ 2.8 L13 | stopper drill | ∅ 4.2 L13 |
| stopper drill | ∅ 3.2 L6 | stopper drill | ∅ 5.2 L6 |
| stopper drill | ∅ 3.2 L8 | stopper drill | ∅ 5.2 L8 |
| stopper drill | ∅ 3.2 L10 | stopper drill | ∅ 5.2 L10 |
| stopper drill | ∅ 3.2 L11.5 | stopper drill | ∅ 5.2 L11.5 |
| stopper drill | ∅ 3.2 L13 | stopper drill | ∅ 5.2 L13 |





Professional Surgical Kit Contain:

Torque ratchet 6.35 hex
Depth probe till 16mm length
Marking drill
2 pilot drill
2.5 drill
2.8 drill
3.2 drill
3.65 drill
4.2 drill
5.2 drill
3.75-4.2 countersink drill
5-6 countersink drill
Drill extender
Guided space pin

BIO-X1021
BIO-X1025
BIO-D3410
BIO-D1220
BIO-D1225
BIO-D1228
BIO-D1232
BIO-D1236
BIO-D1242
BIO-D1252
BIO-D1034
BIO-D1056
BIO-D3142
BIO-GUPI

Implant motor mount 2.42 hex short
Implant motor mount 2.42 hex long
Prosthetics motor mount 1.25 hex
Prosthetics driver 1.25 hex short
Prosthetics driver 1.25 hex long
Implant driver 2.42 hex short
Implant driver 2.42 hex long
Slim implant driver 2.1 hex regular
Conical implant driver regular
X2 parallel pin short
X2 parallel pin long
One-piece implant driver regular
Prosthetic hand driver long

BIO-X1015
BIO-X1014
BIO-X1008
BIO-X1207
BIO-X1215
BIO-X2607
BIO-X2615
BIO-X2810
BIO-X2315
BIO-X1028
BIO-X1029
BIO-OPKEY
BIO-X1006



Pterygoid kit



Pterygoid surgical kit contain:

Short implant driver 2.42mm hex
Medium implant driver 2.42mm hex
Long implant driver 2.42mm hex
Short prosthetic's driver 1.25mm hex
Medium prosthetic's driver 1.25mm hex
Long prosthetic's driver 1.25mm hex
15 degrees dummy abutment
25 degrees dummy abutment
35 degrees dummy abutment
45 degrees dummy abutment
Marking drill sharp
Pterygoid drill 2 step 1.2-2.0
Pterygoid drill 2 step 2.5-2.8
Pterygoid drill 2 step 3.2-3.65
Surgical screw long driver
Pterygoid depth probe 25mm

BIO-X2607
BIO-X2610
BIO-X2622
BIO-X1207
BIO-X1210
BIO-X1215
BIO-A3209
BIO-A3409
BIO-A3509
BIO-A4509
BIO-D3410
BIO-PTDRILL1520
BIO-PTDRILL2528
BIO-PTDRILL3236
BIO-X1023
BIO-PTR01





Zygomatic surgical kit contain:

Marking drill sharp
Zygoma drill 3 step long 2-2.5-3mm
Zygoma drill 3 step long 2-2.5-3mm for 3.75mm implant
Zygoma drill 2 step long 2.8-3.2mm
Zygoma drill long 3.2mm
Zygoma drill 3 step short 2-2.5-3mm
Zygoma drill 3 step short 2-2.5-3mm for 3.75mm implant
Zygoma drill 2 step short 2.8-3.2mm
Zygoma drill short 3.2mm
Diamond bur fine
Diamond bur medium
Diamond bur coarse
Surgical screw driver for 2.42mm hex
Surgical screw driver for 6.35mm hex
Zygoma depth probe 60mm

BIO-D3410
BIO-ZIGDRILL2253/L
BIO-ZIGDRILL2253/L
BIO-ZIGDRILL2832/L
BIO-ZIGDRILL32/L
BIO-ZIGDRILL2253/S
BIO-ZIGDRILL2253/S
BIO-ZIGDRILL2832/S
BIO-ZIGDRILL32/S
BIO-DBF
BIO-DBM
BIO-DBC
BIO-X1022
BIO-X1023
BIO-SPR01



Implant Motor

The new generation of power and control



On Request:
AVAILABLE WITH ATS 20:1 CONTRA ANGLE



- P «Programs»**
Scrolls through the programs
Choice of 5 programs settings
Changed program automatically stored
even after switching off -
- Speed**
Sets the rpm. Motor speeds from 20 to
1500 rpm with 20:1 contra-angle
- Torque**
Changes the torque according to the
reduction ratio
- Flow**
3 different flow rates are offered (Min-
Med-Max (Max.= 100ml/min)

Gear Ratio
Sets the reduction of the contra-angle
Options: 1:1 – 16:1 – 20:1 – 24:1 –
32:1 – 64:1 – 80:1


Calibration
Compensates for loss of torque caused by the
contra-angle

FWD/REV
Switches between forward and
reverse drive
Acoustic signal when reverse drive
is activated

Motor on/off
Activate/Deactivate the Motor

Pump on/off
Activate/Deactivate the Pump

- Brushless Micromotor
- Electronic Calibration of ContraAngle
- Available with 3 different Pedal
- Built to Last
- User Friendly
- Torque: 51 Ncm
- 30000 RPM



omega implant is a complete system, specialty developed under consideration of all the needs in implantology and surgery: powerful performance and precise functionality. see for yourself!
complete, powerful and precise!
powerandaccurate

omega implant is a small babletop unit for implantology that combines the power of a brushless motor with a user friendly. all necessary treatment stages of implantology can be carried out safely and precisely using just one contra-angle 20:1.

omega implants is equipped with contra-angle calibration: the precision of the motor is adapted to the contra-angle connected. omega implant also delivers the higheststandards of ergonomics and flexibility. the motor torque is 50 ncm and the motor speed range is 400 - 30,000 rpm (20-1500 with 20: contra-angle).

Implant Motor with integrated coolant pump easy and uncomplicated to use Brushless motor

AVAILABLE WITH 3 DIFFERENT FOOT CONTROL



«DeLuxe» Multifunction Foot Control

Programs switching - Activation and Deactivation of the coolant pump - Clockwise / Anticlockwise rotation - ON/OFF function. The reostatic effect is ensured by means of a specialized «pressure sensor». There are NO mechanical or moving parts.



«Standard» Multifunction Foot Control

Programs switching - Activation and Deactivation of the coolant pump - Clockwise / Anticlockwise rotation - ON/OFF function. with «ERC» (Electronic Reostatic Effect)



«Easy» Foot Control

ON/OFF button with «ERC» (Electronic Reostatic Effect)



Advanced Irrigation Pump

Easy set up of irrigation tubes - Very ergonomic opening and closing system - Better control of irrigation flow with very constant flow rate - Stable and quiet operations

3 Flow Rates available: MIN - MED - MAX

| Model | Code | Description |
|---------------|--------------|--|
| OMEGA IMPLANT | OMG300E/S/DL | Implantology Unit: - Electronic Control Console - Not Autoclavable Brushless Micromotor - «Easy» - «Standard» - DeLuxe» Foot Control - Solution Holder - Motor Holder - 1 Tube |









Open System

output STL file, compatible with 3rd party software. portable design with suitcase. customized PTY format for both clinics and labs.

Application

-  orthodontics
-  implantation
-  prosthodontics
-  cosmetic



High Speed



Easy to handle and Operate



True Color



Applicable for implant& orthodontic



Powder-free



High Accuracy



Replaceable Header



Automatic Preparation Evaluation

Parameters

| | | | |
|---------------------------------|------------------|-----------------------|----------------------------|
| scanheader (L*W*H) ¹ | 110 mm*19mm*17mm | interface | US3.0 |
| size (L*W*H) | 221mm*54mm*17mm | weight | 350 g |
| scanning speed | <3 min | data format | stl, ply, pty ² |
| scanning depth | 0~15mm | accuracy | 15 μm |
| window size | 17mmx20mm | power | adaptor ac 100v~240v |
| cable length | 1.5m | operating environment | 0°C~50°C rh<80% |
| | | storage environment | -10°C~65°C rh<90% |



Digital components



one digital abutment + t base + scan body

| | | | |
|---------|-----------|-----------|-----------|
| Cat No. | BIO-IS-05 | BIO-IS-15 | BIO-IS-25 |
| Height | 0.5 mm | 1.5 mm | 2.5 mm |



t base with hex

| | |
|---------|---------|
| Cat No. | BIO-TBH |
| Height | 5.0 mm |



t base non hex

| | |
|---------|---------|
| Cat No. | BIO-TBN |
| Height | 5.0 mm |



conical t base

| | |
|---------|---------|
| Cat No. | BIO-CTB |
| Height | 5.0 mm |



Bioline Dental GmbH &Co.KG

Akazien str 7 16356 Wermeuchen Germany
Tel +49 15202643207
info@bioline-implants.de

EC **REP** Obelis ^{SA}

European Authorised Representative Center Registered
Adress:

Bd. General Wahis 53 1030 Brussels-Belgium
www.obelis.net
Tel: +32 (0)2 732 59 54

India Head Office:
BioLine Dental Pvt Ltd.
Office Number 27 1st Floor
Shreeji Villa CHS Ltd
T.H Kataria Marg
Matung West
Mumbai 400 016
India.
Tel:+91-9137872067

Customer Support:
info@bioline-implants.de
www.bioline-implants.de

WWW.BIOLINE-IMPLANTS.DE
